

PROF. MARIANNE VAN HAGE (Orcid ID : 0000-0003-3091-1596) DR. HANS GRÖNLUND (Orcid ID : 0000-0003-4882-7624) PROF. RUDOLF VALENTA (Orcid ID : 0000-0001-5944-3365) DR. MIRELA CURIN (Orcid ID : 0000-0003-2575-4127)

Article type : Letter to the Editor

# Highly sensitive ELISA-based assay for quantification of allergen-specific IgE antibody levels

Antonina Karsonova<sup>1</sup>, Ksenja Riabova<sup>1</sup>, Sergio Villazala-Merino<sup>2</sup>, Raffaela Campana<sup>3</sup>, Verena Niederberger<sup>2</sup>, Julia Eckl-Dorna<sup>2</sup>, Renate Fröschl<sup>4</sup>, Thomas Perkmann<sup>4</sup>, Yury V. Zhernov<sup>5</sup>, Olga G. Elisyutina<sup>5</sup>, Elena S. Fedenko<sup>5</sup>, Musa R. Khaitov<sup>5</sup>, Daria Fomina<sup>1,6</sup>, Evgeniy Beltiukov<sup>7</sup>, Marianne van Hage<sup>8</sup>, Hans Grönlund<sup>9</sup>, Rudolf Valenta<sup>3,10</sup>, Alexander Karaulov<sup>1</sup>, Mirela Curin<sup>3</sup>

<sup>1</sup>Laboratory of Immunopathology, Department of Clinical Immunology and Allergy, Sechenov First Moscow State Medical University, Moscow, Russian Federation. <sup>2</sup>Department of Otorhinolaryngology, Medical University of Vienna, Vienna, Austria. <sup>3</sup>Division of Immunopathology, Department of Pathophysiology and Allergy Research, Center for Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Vienna, Austria <sup>4</sup>Department of Laboratory Medicine, Medical University of Vienna, Vienna, Austria. <sup>5</sup>National Research Center - Institute of Immunology FMBA of Russia, Moscow, Russian Federation. <sup>6</sup>City Hospital #52, Moscow, Russian Federation. <sup>7</sup>Ural State Medical University, Ekaterinburg, Russian Federation. <sup>8</sup>Division of Immunology and Allergy, Department of Medicine Solna, Karolinska Institutet and University Hospital, SE-171 77, Stockholm, Sweden. <sup>9</sup>Department of Clinical

This article has been accepted for publication and undergone full peer review but has not been through the copyediting, typesetting, pagination and proofreading process, which may lead to differences between this version and the <u>Version of Record</u>. Please cite this article as <u>doi:</u> 10.1111/ALL.14325

This article is protected by copyright. All rights reserved

Neuroscience, Therapeutic Immune Design Unit, Karolinska Institutet, Stockholm, Sweden. <sup>10</sup>Karl Landsteiner University for Healthcare Sciences, Krems, Austria.

#### **Correspondence:**

Mirela Curin, Department of Pathophysiology and Allergy Research, Center for Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Währinger Gürtel 18-20, A-1090 Vienna, Austria. Email: mirela.curin@meduniwien.ac.at

Key words: Allergy, allergen, IgE quantification, recombinant allergen, standard curve.

### To the Editor:

Allergy is the most common immunologically-mediated hypersensitivity disease worldwide affecting upper and lower respiratory tract, the skin as well as the gastrointestinal tract and may cause life-threatening anaphylaxis<sup>1</sup>. Immunoglobulin E (IgE) occurs in extremely low amounts in blood and tissues but is the key immunoglobulin in allergy. The quantification of allergen-specific IgE levels (sIgE) in serum and body fluids is important because sIgE concentrations are useful to define clinical phenotypes of allergy and to predict the intensity of allergic reactions<sup>2-5</sup>. However, quantification of sIgE requires special instruments (e.g., ImmunoCAP) that are available only in industrialized countries, large volumes of serum and is quite expensive<sup>6</sup>. Therefore, the goal of our study was to establish a quantitative method for measuring allergen-specific IgE based on a simple ELISA technology which can be implemented and used in all parts of the world. The user will thus need only equipment for ELISA, reagents, detection antibodies and test allergens plus reagents for the IgE standard curve.

In order to establish a standard for ELISA-based quantification of allergen-specific human IgE, a human monoclonal chimeric IgE antibody (IgEmoAb)<sup>7</sup> consisting of the variable region of a mouse  $IgG_1$  antibody specific for the major birch pollen allergen Bet v 1 which had been fused to the human epsilon heavy chain was expressed and purified from hybridoma cells by affinity chromatography to >95% purity (Figure S 1A). Under non-reducing conditions purified IgEmoAb

migrates as a distinct single band of approximately 250 kDa molecular weight in SDS-PAGE, while under reducing conditions two distinct bands which correspond to the heavy and light chain became visible at approximately 100 kDa and 20 kDa, respectively (Figure S 1A). IgEmoAb reacted specifically with anti-human IgE antibodies under non-reducing conditions and under reducing conditions the heavy chain showed selective anti-IgE reactivity by immunoblotting (Figure S 1B). Structural analysis of IgEmoAb by circular dichroism showed that it assumed a predominant beta sheet fold with a minimum signal around 217 nm and a maximum near 198 nm (Figure S 1C).

Using IgEmoAb and Bet v 1 we next established by ELISA a standard curve allowing to measure and quantify sIgE. A series of dilutions of IgEmoAb ranging from 2.42ng/ml to 2420ng/ml was exposed to ELISA plate-bound IgE yielding standard curves with a linear range of approximately 2.42 – 242 ng/ml IgE depending on the time of assay development (Figure 1A). Supplementary Figure 2 shows the OD values obtained for the dilution series of IgEmoAb determined by ELISA at different time points of assay development i.e., (after 5 min, 10 min, 15 min and 20 min) yielding optimal results from 5-10 minutes of development with the given detection antibodies and color reagents. The different dilutions of IgEmoAb were also measured using Streptavidin ImmunoCAPs prepared with the same recombinant Bet v 1 preparation used in the ELISA experiments. The concentrations of IgEmoAb i.e., 2.42 ng/ml, 24.2 ng/ml, 242 ng/ml, 2420 ng/ml corresponded to the following ImmunoCAP results: 0.34 kU<sub>A</sub>/L, 3.28 kU<sub>A</sub>/L, 33.1 kU<sub>A</sub>/L and  $>100 \text{ kU}_{\text{A}}/\text{L}$  showing that the dilutions of the IgEmoAb standard were correct down to few nanograms per ml without loss of activity. However, we noted that the ImmunoCAP determination of the concentrations of IgEmoAb were different (i.e., yielding approx. 30%) from the results obtained by microBCA protein determination which may be due to differences between the ImmunoCAP IgE standard and IgEmoAB (data not shown). Accordingly a "real conversion factor" of 7ng/ml = 1 kUA/L for calculating IgE concentrations had to be introduced. Based on a standard curve established with IgEmoAb, Bet v 1-specific IgE concentrations were quantified in 19 sera from clinically well-characterized birch pollen allergic patients by ELISA (ng/ml) (Tables S1 and S2). Table S2 displays the OD levels, the corresponding ng/ml concentrations of Bet v 1-specific IgE and the results from the ImmunoCAP and ImmunoCAP ISAC measurements in detail. We noted a strong significant correlation between the Bet v 1specific IgE concentrations determined in ng/ml with the ELISA and in kU<sub>A</sub>/l ImmunoCAP results (r=0.904, p<0.0001, 95% CI 0.756 to 0.964; Figure 1B). In fact, Bet v 1-specific IgE levels

measured by ELISA correlated even better with ImmunoCAP results than IgE levels measured with ImmunoCAP ISAC chip technology (r=0.875, p<0.0001, 95% CI 0.692 to 0.953; Figure 1C). Nevertheless, results of the quantitative IgE ELISA correlated also significantly with sIgE levels measured with ImmunoCAP ISAC chip technology (r=0.840, p<0.0001, 95% CI 0.616 to 0.939) (Figure 1D). In summary, strong correlations of Bet v 1-specific IgE levels measured by ELISA with those determined by ImmunoCAP and ImmunoCAP ISAC technology were obtained.

Next, we investigated whether the IgEmoAb-based standard curve is also suitable for the quantification of patients' IgE to an allergen molecule unrelated to Bet v 1. For this purpose we selected the major cat allergen Fel d 1 which besides Bet v 1 is one of the most frequently recognized allergens in the Russian population<sup>5,8</sup>. Quantitative IgE ELISA was used to determine Fel d 1-specific IgE levels in 19 cat allergic patients using different serum dilutions. Detailed results of quantitative IgE ELISA and of ImmunoCAP are shown in Table S3 and Table S4. A strong significant correlation was obtained between Fel d 1-specific IgE concentrations quantified with the IgE ELISA (ng/ml) and with the ImmunoCAP results(kU<sub>A</sub>/l) (r=0.984, p<0.0001, 95% CI 0.956 to 0.994; Figure 2A) demonstrating that the quantitative ELISA can be used for quantification of specific IgE levels against different allergens, e.g., Fel d 1. Interestingly, we noted for three patients with clinically confirmed symptoms of cat allergy ImmunoCAP measurements were below the threshold of 0.1 kUA/L but we were able to detect and quantify Fel d 1-specific IgE with the IgE ELISA (Table S3, patients C5, C6 and C13). Of note, we found that the results obtained for different dilutions of the sera delivered highly comparable results. Results of IgE level assessment obtained from 1:2 sera dilutions and 1:4 sera dilutions from the cat allergic patients (shown in Table S4) were compared and a strong and significant correlation was noted (r= 0.988 p<0.0001, 95% CI 0.967 to 0.996; Figure 2B).

The quantitative IgE ELISA assay should be useful for laboratories in less developed countries which cannot afford expensive IgE test systems for routine allergy diagnosis taking advantage of molecular allergy diagnosis<sup>9,10</sup>. Before extensive use, broader clinical studies correlating IgE concentrations to clinical signs may be warranted. In fact, the price for a quantitative IgE determination for one allergen should cost less than 10% than by ImmunoCAP, however, the time required for testing will be longer. Our assay will therefore be useful in places where work power cost is not a major issue and the results are desired in a higher-throughput format as well as for cases (e.g., research) where regulations are of minor importance. Furthermore, the assay should

allow a flexible testing of different allergen molecules which can be kept on stock and used depending on the clinical questions and diagnostic needs required. To achieve these goals the authors will make efforts to make available the key reagents for the standard curve (i.e., IgEmoAb and rBet v 1) to the community. In summary, we provide a simple, inexpensive, sensitive and accurate method for measuring allergen-specific IgE concentrations in serum and body fluids to facilitate the diagnosis of IgE-associated allergy in the world.

Antonina Karsonova<sup>1</sup> Ksenja Riabova<sup>1</sup> Sergio Villazala-Merino<sup>2</sup> Raffaela Campana<sup>3</sup> Verena Niederberger<sup>2</sup> Julia Eckl-Dorna<sup>2</sup> Renate Fröschl<sup>4</sup> Thomas Perkmann<sup>4</sup> Yury V. Zhernov<sup>5</sup> Olga G. Elisyutina<sup>5</sup> Elena S. Fedenko<sup>5</sup> Musa R. Khaitov<sup>5</sup> Daria Fomina<sup>1,6</sup> Evgeniy Beltiukov<sup>7</sup> Marianne van Hage<sup>8</sup> Hans Grönlund<sup>9</sup>

Rudolf Valenta<sup>3, 10</sup>

Alexander Karaulov<sup>1</sup>

Mirela Curin<sup>3</sup>

# From:

<sup>1</sup>Laboratory of Immunopathology, Department of Clinical Immunology and Allergy, Sechenov
First Moscow State Medical University, Moscow, Russian Federation. <sup>2</sup>Department of
Otorhinolaryngology, Medical University of Vienna, Vienna, Austria. <sup>3</sup>Division of
Immunopathology, Department of Pathophysiology and Allergy Research, Center for
Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Vienna, Austria
<sup>4</sup>Department of Laboratory Medicine, Medical University of Vienna, Vienna, Austria. <sup>5</sup>National
Research Center - Institute of Immunology FMBA of Russia, Moscow, Russian Federation. <sup>6</sup>City
Hospital #52, Moscow, Russian Federation. <sup>7</sup>Ural State Medical University, Ekaterinburg, Russian
Federation. <sup>8</sup>Division of Immunology and Allergy, Department of Medicine Solna, Karolinska
Institutet and University Hospital, SE-171 77, Stockholm, Sweden. <sup>9</sup>Department of Clinical
Neuroscience, Therapeutic Immune Design Unit, Karolinska Institutet, Stockholm, Sweden. <sup>10</sup>Karl
Landsteiner University for Healthcare Sciences, Krems, Austria.

# **References:**

1.

2.

- Anto, J. M. *et al.* Mechanisms of the Development of Allergy (MeDALL): Introducing novel concepts in allergy phenotypes. *J. Allergy Clin. Immunol.* **139**, 388–399 (2017).
- Westman, M. *et al.* Early childhood IgE reactivity to pathogenesis-related class 10 proteins predicts allergic rhinitis in adolescence. *J. Allergy Clin. Immunol.* **135**, 1199-206.e1–11 (2015).
  - Wickman, M. *et al.* Detection of IgE reactivity to a handful of allergen molecules in early childhood predicts respiratory allergy in adolescence. *EBioMedicine* **26**, 91–99 (2017).
- 4. Grönlund, H. *et al.* Higher immunoglobulin E antibody levels to recombinant Fel d 1 in catallergic children with asthma compared with rhinoconjunctivitis. *Clin. Exp. Allergy* **38**,

1275-81 (2008).

- Elisyutina, O. *et al.* Bet v 1-specific IgE levels and PR-10 reactivity discriminate silent sensitization from phenotypes of birch allergy. *Allergy* **74**, 2525–2528 (2019).
- van Hage, M., Hamsten, C. & Valenta, R. ImmunoCAP assays: Pros and cons in allergology. *J. Allergy Clin. Immunol.* **140**, 974–977 (2017).
- Laffer, S. et al. A molecular model of type I allergy: identification and characterization of a nonanaphylactic anti-human IgE antibody fragment that blocks the IgE-FcepsilonRI interaction and reacts with receptor-bound IgE. J. Allergy Clin. Immunol. 108, 409–16 (2001).
- Grönlund, H. *et al.* Formation of disulfide bonds and homodimers of the major cat allergen Fel d 1 equivalent to the natural allergen by expression in Escherichia coli. *J. Biol. Chem.* **278**, 40144–51 (2003).
- Levin, M. E. *et al.* Environmental factors associated with allergy in urban and rural children from the South AFrican Food Allergy (SAFFA) cohort. *J. Allergy Clin. Immunol.* (2019). doi:10.1016/j.jaci.2019.07.048
- Gray, C. L., Levin, M. E. & du Toit, G. Egg sensitization, allergy and component patterns in African children with atopic dermatitis. *Pediatr. Allergy Immunol.* **27**, 709–715 (2016).

Acknowledgments: This study was supported by a Megagrant of the Government of the Russian Federation, grant number 14.W03.31.0024, by the "Russian Academic Excellence Project 5-100" and by the Austrian Science Fund (FWF) projects F4605, DK1248-B13, and F4613.

# **Author Contributions**

R.V. designed, supervised and coordinated the study. M.C, AN.K, K.R. and R.V. participated in the interpretation of the findings and wrote the manuscript. AN.K., K.R., S.V.M., J.E.D., R.C., R.F., Y.Z. and O.E., conducted experiments, analyzed the data and prepared the figures. M.K.,

AL.K., V.N., T.P., D.F., E.B., P.G., M.V.H., H.G. and E.F. recruited and characterized the patients, collected the sera, analyzed and interpreted the results. All authors revised and approved the manuscript.

### **Competing interests**

R.V. has received research grants from HVD company, Biomay AG, Vienna, Austria and Viravaxx, Vienna, Austria and serves as consultant for Viravaxx. MvH has received lecture fees from Thermo Fisher Scientific and consultancy fees from Biomay AG, Vienna Austria and Hycor Biomedical LLC, CA, US outside of the submitted work. Other authors have no conflict of interest to declare.

#### **Figure legends**

**Figure 1.** Association of different concentrations of IgEmoAb (x-axis: ng/ml) with IgE reactivities to Bet v 1 determined as optical density (OD) values (y-axis) by ELISA (**A**). Correlation of Bet v 1-specific IgE levels determined in sera from 19 birch pollen allergic patients by different methods. (**B**) ImmunoCAP ( $kU_A/L$ ) (y-axis) versus ng/ml determined by ELISA (x-axis). (**C**) ISU (ISAC standardized unit) values determined by ImmunoCAP ISAC (y-axis) versus ImmunoCAP ( $kU_A/L$ ) (x-axis) and (**D**) ISU values determined by ImmunoCAP ISAC (y-axis) versus ng/ml determined by ELISA (x-axis). Correlation coefficients (r), regression equations, and statistical significances (p values) are shown in the top upper corner of graph.

**Figure 2**. Correlation of Fel d 1-specific IgE levels determined in sera from 19 cat allergic patients by ImmunoCAP ( $kU_A/L$ ) (y-axis) versus ng/ml determined by ELISA (x-axis) (**A**). Correlation of Fel d 1-specific IgE levels determined in lower (y-axis) and higher (x-axis) sera dilutions from 19 cat allergic patients from Table S3 (**B**).Correlation coefficients (r), regression equations, and statistical significances (p values) are shown in the top upper corner of graph.





Figure 2